



QUANTA BIODESIGN LIMITED

COO: US 7470 Montgomery Drive, Plain City, OH 43064
Phone: (614) 792-2958 Fax: (614) 760-9781 www.quantabiodesign.com

62-0200-0800

Q-Bright Endure Plus Chemiluminescent Detection Kit 62-0200-0800

Q-Bright Endure Plus (62-0200-0800) chemiluminescent detection reagent is a luminol-based substrate for horseradish peroxidase (HRP), providing a light output for up to 24h. *Q-Bright* Endure Plus provides 2X the signal of *Q-Bright* Endure without compromising on signal duration, providing long-lasting, mid-femtogram levels of detection in western blot, ELISA and other chemiluminescent-based applications.

Components:	32-0208-0200 Chemiluminescent Substrate Solution 32-1221-0200 Stable Peroxide Solution
Sensitivity:	mid-femtogram
Preservative:	0.1% Kathon CG/ICP
Storage and Handling:	Store at 4° C. <i>May ship at ambient temperature.</i>
Unit sizes:	<i>Q-Bright</i> Endure Plus is available in 20ml, 100ml, and 200ml unit sizes.

Instructions for Use

Q-Bright Endure Plus Chemiluminescent Detection Kit (62-0200-0800) Working Solution is prepared by combining equal parts of the Chemiluminescent Substrate Solution and the Stable Peroxide Solution before use. Allow working solution to equilibrate to room temperature before applying to the membrane.

For western blots, using optimized concentrations of primary antibody (0.02-0.5ug/ml) and HRP conjugated secondary antibody (2-20ng/ml), develop the membrane with a minimum of 0.1ml of Working Solution per cm² of membrane. Allow signal to develop for 5 minutes before imaging. The chemiluminescent signal can be developed using x-ray film or chemiluminescent imagers.

Chemiluminescent Product Line

The *Q-Bright* series of chemiluminescent reagents for HRP detection come in a variety of formulations and working ranges. When matched with the appropriate Antibody-HRP or Streptavidin-HRP conjugate, the result is high sensitivity, low background and sustained duration in immunoassays.