



Quanta BioDesign, Ltd

7500 Montgomery Drive,
Plain City, Ohio 43064

Office 614.792.2958

Fax 614.760.9781

QuantaBioDesign.com

AQ LIFE SCIENCES
A DIVISION OF QUANTA BIODESIGN LIMITED

Product Number: AQ-30024

Goat anti-Mouse IgG (H&L), dPEG[®] Horseradish Peroxidase min X Rabbit, Human, Bovine, Horse *min X Rabbit, Human, Bovine, Horse*

dPEG[®] Horseradish Peroxidase (HRP) conjugated Goat anti-Mouse IgG (H&L) adsorbed to rabbit, human, bovine and horse IgG is useful in immunodetection applications. Goat anti-Mouse IgG (H&L) has been covalently conjugated to horseradish peroxidase (HRP) utilizing our exclusive SuperHydrophilic technology designed to improve solubility and limit non-specific binding. HRP-conjugated Goat anti-Mouse IgG (H&L) with this adsorption profile is useful in the classic range of applications including immunohistochemistry, western blots, and ELISAs using multiple labels. For best results, use this Goat anti-Mouse IgG (H&L), dPEG[®] HRP conjugate in conjunction with our Q-Bright[®] Chemiluminescence detection reagents for ELISA and western blotting applications.

Product Number:

AQ-30024

Concentration:

1 mg/ml

Buffer:

PBS, pH 7.2, with 1% BSA, and other proprietary stabilizers

Preservative:

0.1% Kathon CG/ICP

Storage and Handling:

Store at 4°C for up to 12 months. To extend shelf life, aliquot and store at -20°C or below and avoid repeated freeze-thaw cycles. Alternatively, add an equal volume of molecular biology grade glycerol (50% final concentration) and store at -20°C as a liquid. Protect from light.

Antibody Information

Target Antigen:

Mouse IgG (H&L)

Host:

Goat

Clonality:

Polyclonal

Minimal Cross Reactivity to:

Rabbit, Human, Bovine, Horse

Antibody Country of Origin:

United States of America

Label Information

dPEG[®] HRP

FOR RESEARCH USE ONLY. THIS PRODUCT IS A LABORATORY REAGENT AND IS NOT TO BE ADMINISTERED TO HUMANS OR ANIMALS. IT IS NOT TO BE USED FOR ANY DIAGNOSTIC OR THERAPEUTIC PURPOSE. NOT FOR RESALE.

SuperHydrophilic™ is a trademark of Quanta BioDesign, Ltd. dPEG[®], AQuora[®], Q-Bright[®], AQ[®], Q-Link[®], and Q-Tag[®] are all registered trademarks of Quanta BioDesign, Ltd.